

Product Datasheet

phospho-p38 MAPK (Thr180 + Tyr182) Antibody (BF750) (orb1593745)

Description

phospho-p38 MAPK (Thr180 +

Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	BF750
Tested Applications	FC, ICC
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human p38 MAPK around the phosphorylation site of Thr180/Tyr182
Target	phospho-p38 MAPK (Thr180 + Tyr182)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q16539
Dilution Range	Flow-Cyt=1µg /test,ICC=1:100
Expiration Date	12 months from date of receipt.