

Product Datasheet

phospho-IRS1 (Ser636 + Ser639) Antibody (PerCP) (orb1594314)

Description	phospho-IRS1 (Ser636 + Ser639) Antibody (PerCP)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	PerCP
Tested Applications	FC, ICC
Immunogen	KLH conjugated synthesised phosphopeptide derived from human IRS1 around the phosphorylation site of Ser636/639
Target	phospho-IRS1 (Ser636 + Ser639)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P35568
Dilution Range	Flow-Cyt=1µg /test,ICC=1:100
Expiration Date	12 months from date of receipt.