

Product Datasheet

Phospho-P38 MAPK (Thr180 + Tyr182) Antibody (BF680) (orb1594697)

Description	Phospho-P38 MAPK (Thr180 + Tyr182) Antibody (BF680)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	BF680
Tested Applications	FC, IF
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human p38 MAPK around the phosphorylation site of Thr180/Tyr182
Target	Phospho-P38 MAPK (Thr180 + Tyr182)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q16539
Dilution Range	Flow-Cyt=1µg/Test,IF=1:100-500
Expiration Date	12 months from date of receipt.