

Product Datasheet

JNK1+JNK2+JNK3 Antibody (BF405) (orb1594863)

Description	JNK1+JNK2+JNK3 Antibody (BF405)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	BF405
Tested Applications	FC, ICC, IF
Immunogen	KLH conjugated synthetic peptide derived from human JNK1/2/3
Target	JNK1+JNK2+JNK3
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P45983
Dilution Range	Flow-Cyt=1ug/Test,IF=1:100-500,ICC=1:100
Expiration Date	12 months from date of receipt.