

Product Datasheet

KCNQ2 Antibody (Cy5.5) (orb1595095)

Description	KCNQ2 Antibody (Cy5.5)
Species/Host	Rabbit
Reactivity	Rat
Conjugation	PerCP/Cy5.5
Tested Applications	FC, ICC
Immunogen	KLH conjugated synthetic peptide derived from human KCNQ2
Target	KCNQ2
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	043526
Dilution Range	Flow-Cyt=3µg/Test,ICC=1:100-500
Expiration Date	12 months from date of receipt.