

Product Datasheet

GPCR ORF4/GPCR GPR137 Antibody (BF680) (orb1595584)

Description	GPCR ORF4/GPCR GPR137 Antibody (BF680)
Species/Host	Rabbit
Reactivity	Human
Conjugation	BF680
Tested Applications	ICC, IF
Immunogen	KLH conjugated synthetic peptide derived from human GPCR ORF4/GPCR GPR137
Target	GPCR ORF4/GPCR GPR137
Preservatives	Preservative: 15mM Sodium Azide, Constituents: 1% rAlbumin, 0.01M PBS, pH 7.4
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q96N19
Dilution Range	IF=1:100-500,ICC=1:100-500
Expiration Date	12 months from date of receipt.