

Product Datasheet

Thyroid hormone receptor alpha Antibody (AP) (orb1596346)

Description	Thyroid hormone receptor alpha Antibody (AP)
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	AP
Tested Applications	ELISA, WB
Immunogen	KLH conjugated synthetic peptide derived from human Thyroid hormone receptor alpha
Target	Thyroid hormone receptor alpha
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P10827
Dilution Range	WB=1:500-2000,ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.