

Product Datasheet

Vinculin (Loading Control) Antibody (Biotin) [Out of stock] (orb1596500)

Description	Vinculin (Loading Control) Antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Biotin
Tested Applications	WB
Immunogen	Recombinant Human Vinculin
Target	Vinculin (Loading Control)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Uniprot ID	P18206
Dilution Range	WB: 1:500-2000, ICC: 1:50, IF: 1:50-200
Expiration Date	12 months from date of receipt.