

Product Datasheet

Beta Actin (Yeast, Loading Control) Antibody (AP) [Out of stock] (orb1597014)

Description	Beta Actin (Yeast, Loading Control) Antibody (AP)
Species/Host	Rabbit
Reactivity	Yeast
Conjugation	AP
Tested Applications	WB
Immunogen	Recombinant Yeast Actin.
Target	Beta Actin (Yeast, Loading Control)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P60010
Dilution Range	WB: 1:1000-5000
Expiration Date	12 months from date of receipt.