

Product Datasheet

NF-L Antibody (Biotin) (orb1597143)

Description	NF-L Antibody (Biotin)
Species/Host	Mouse
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA, WB
Immunogen	Recombinant Human NF-L Protein
Target	NF-L
Preservatives	0.01M PBS(pH7.4)
Form/Appearance	Liquid
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Uniprot ID	P07196
Dilution Range	WB=1:500-2000,ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.