


Product Datasheet

CD52 Antibody (orb1597172)

Description	CD52 Antibody
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human CD52 (20-50/61 aa)
Target	CD52
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	6.6 kDa
Uniprot ID	P31358
Dilution Range	IHC-P=1:50-200
Expiration Date	12 months from date of receipt.



IHC-P
analysis of
rat spleen
using CD52

...