

Product Datasheet

CH25H Antibody (HRP) (orb1597314)

Description	CH25H Antibody (HRP)
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	HRP
Tested Applications	WB
Immunogen	Recombinant human CH25H protein
Target	CH25H
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Lyophilized or Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O95992
Dilution Range	IHC-P=1:100-500,ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.