

Product Datasheet

PLK1 Antibody (Biotin) (orb1597593)

Description	PLK1 Antibody (Biotin)
Species/Host	Mouse
Reactivity	Human
Conjugation	Biotin
Tested Applications	WB
Immunogen	Recombinant human PLK1 protein
Target	PLK1
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
CAS Number	Serine/threonine-protein kinase PLK1
Uniprot ID	P53350
Expiration Date	12 months from date of receipt.