

Product Datasheet

PPV-VP2 Antibody (orb1597662)

Description	PPV-VP2 Antibody
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant PPV VP2 protein (1-579/579 aa)
Target	PPV-VP2
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
MW	64 kDa
Uniprot ID	P18546
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.

