

Product Datasheet

GPR48 Antibody (PerCP) (orb1597915)

Description	GPR48 Antibody (PerCP)
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	PerCP
Tested Applications	FC
Immunogen	KLH conjugated synthetic peptide derived from human LGR4
Target	GPR48
Preservatives	Preservative: 15mM Sodium Azide, Constituents: 1% rAlbumin, 0.01M PBS, pH 7.4
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9BXB1
Expiration Date	12 months from date of receipt.