

Product Datasheet

SARS-CoV-2 N Antibody (HRP) (orb1599854)

Description	SARS-CoV-2 N Antibody (HRP)
Species/Host	Rabbit
Reactivity	Virus
Conjugation	HRP
Tested Applications	WB
Immunogen	Recombinant SARS-CoV-2 Nucleocapsid protein
Target	SARS-CoV-2 N
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Lyophilized or Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	46 kDa
Uniprot ID	P0DTC2
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.