

Product Datasheet

SARS-CoV Nucleocapsid Antibody (HRP) (orb1599863)

Description	SARS-CoV Nucleocapsid Antibody (HRP)
Species/Host	Mouse
Conjugation	HRP
Tested Applications	ELISA, WB
Immunogen	Recombinant SARS N (nucleocapsid)
Target	SARS-CoV Nucleocapsid
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	≥1mg/ml
Storage	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
Uniprot ID	P59594
Dilution Range	WB=1:500-2000,ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.