


Product Datasheet

Alpha casein Antibody (orb1600109)

Description	Alpha casein Antibody
Species/Host	Rabbit
Reactivity	Bovine
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant bov Alpha s1 Casein protein (16-214/214 aa)
Target	Alpha casein
Preservatives	0.01M TBS (pH7.4) with 1% rAlbumin, 0.02% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	23 kDa
Uniprot ID	P47710
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.



Western blot: Lane 1: Recombinant Bovine A...