

Product Datasheet

VWF Antibody (BF350) (orb1601215)

Description	VWF Antibody (BF350)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	BF350
Tested Applications	FC
Immunogen	KLH conjugated synthetic peptide derived from human VWF
Target	VWF
Preservatives	Preservative: 15mM Sodium Azide, Constituents: 1% rAlbumin, 0.01M PBS, pH 7.4
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P04275
Expiration Date	12 months from date of receipt.