

Product Datasheet

phospho-STAT3 (Tyr705) Antibody (Biotin) (orb1602002)

Description	phospho-STAT3 (Tyr705) Antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Biotin
Tested Applications	WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human STAT3 around the phosphorylation site of Tyr705
Target	phospho-STAT3 (Tyr705)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P40763
Dilution Range	WB=1:200-1000,IHC-P=1:100-500,IHC-F=1:100-500,ICC=1:100-500
Expiration Date	12 months from date of receipt.