


Product Datasheet

ATG5/APG5L Antibody (orb1602288)

Description	ATG5/APG5L Antibody
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC-Fr, WB
Immunogen	Recombinant human ATG5 Protein
Target	ATG5/APG5L
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
MW	32 kDa
Uniprot ID	Q9H1Y0
Dilution Range	WB=1:500-1000, IHC-F=1:100-500, IF=1:100-500, IHC-P=1:100-500
Expiration Date	12 months from date of receipt.



Western blot: Lane 1: Human PANC-1 cell; Lane 2: Human MCF-7 cell