

Product Datasheet

Phospho-TrkA (Tyr496) +TrkB(Tyr516) Antibody (Cy5.5) (orb1602976)

Description	Phospho-TrkA (Tyr496) +TrkB(Tyr516) Antibody (Cy5.5)
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	PerCP/Cy5.5
Tested Applications	FC
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human TrkA around the phosphorylation site of Tyr496
Target	Phospho-TrkA (Tyr496) +TrkB(Tyr516)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P04629
Expiration Date	12 months from date of receipt.