

Product Datasheet

Hepatitis C Virus RNA directed RNA polymerase antibody (orb1604)

Description Hepatitis C Virus RNA directed RNA polymerase antibody

Species/Host Rabbit

Conjugation Unconjugated

Tested Applications ELISA

Immunogen KLH conjugated synthetic peptide derived from Hepatitis C Virus RNA-directed RNA polymerase (2501-2600/3010 aa)

Target Hepatitis C Virus RNA-directed RNA polymerase

Preservatives 0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.

Form/Appearance Liquid

Concentration 1mg/ml

Storage Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.

Note For research use only


Isotype IgG

Clonality Polyclonal

MW 65 kDa

Uniprot ID [Q14SZ0](#)

Expiration Date 12 months from date of receipt.



Western blot analysis of rat intestine I...