

Product Datasheet

GRP78/Bip Antibody (BF594) (orb1604553)

Description	GRP78/Bip Antibody (BF594)
Species/Host	Mouse
Reactivity	Human
Conjugation	BF594
Tested Applications	ICC, IF
Immunogen	KLH conjugated synthetic peptide derived from human GRP78/Bip
Target	GRP78/Bip
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Uniprot ID	P11021
Dilution Range	IF=1:100-500,ICC=1:100
Expiration Date	12 months from date of receipt.