

Product Datasheet

SMAD7 Antibody (HRP) (orb1605230)

Description	SMAD7 Antibody (HRP)
Species/Host	Mouse
Reactivity	Human
Conjugation	HRP
Tested Applications	IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human SMAD7
Target	SMAD7
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Uniprot ID	O15105
Dilution Range	WB=1:500-2000,IHC-P=1:100-500
Expiration Date	12 months from date of receipt.