

Product Datasheet

CK14+17+42+10 Antibody (Biotin) (orb1605334)

Description	CK14+17+42+10 Antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Biotin
Tested Applications	IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human CK14+17+42+10
Target	CK14+17+42+10
Preservatives	0.01M TBS (pH7.4)with 0.05% rAlbumin, 40% Glycerol, 0.05% Sodium Azide.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P08779
Dilution Range	WB=1:500-2000,IHC-P=1:100-500
Expiration Date	12 months from date of receipt.