

Product Datasheet

phospho-MAPK8 (Thr183 + Tyr185) Antibody (Cy5.5) (orb1606396)

Description	phospho-MAPK8 (Thr183 + Tyr185) Antibody (Cy5.5)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	APC/Cy5.5
Tested Applications	FC
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human MAPK9 around the phosphorylation site of Thr183+Tyr185
Target	phospho-MAPK8 (Thr183 + Tyr185)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P45983
Expiration Date	12 months from date of receipt.