

Product Datasheet

phospho-Tau (Thr181) Antibody (BF680) (orb1606595)

Description	phospho-Tau (Thr181) Antibody (BF680)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	BF680
Tested Applications	FC, ICC, IF
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human Tau around the phosphorylation site of Thr181
Target	phospho-Tau (Thr181)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P10636
Dilution Range	Flow-Cyt=1ug/test,IF=1:50-200,ICC=1:100
Expiration Date	12 months from date of receipt.