

Product Datasheet

Phospho-ATG1 (Ser556) Antibody (BF680) (orb1606955)

Description	Phospho-ATG1 (Ser556) Antibody (BF680)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	BF680
Tested Applications	FC, ICC, IF
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human ATG1/ULK1 around the phosphorylation site of Ser556
Target	Phospho-ATG1 (Ser556)
Preservatives	Preservative: 15mM Sodium Azide, Constituents: 1% rAlbumin, 0.01M PBS, pH 7.4
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	075385
Dilution Range	Flow-Cyt=3ug/test,IF=1:100-500,ICC=1:100
Expiration Date	12 months from date of receipt.