

Product Datasheet

pRUNX3 antibody (orb1607)

Description	pRUNX3 antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	KLH conjugated synthesised phosphopeptide derived from human RUNX3 around the phosphorylation site of Ser338 (TS(p-S)GS)
Target	phospho-RUNX3 (Ser338)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	44 kDa
Uniprot ID	Q13761
Expiration Date	12 months from date of receipt.