

Product Datasheet

phospho-GSK-3 Beta (Ser9) Antibody (HRP) (orb1607151)

Description	phospho-GSK-3 Beta (Ser9) Antibody (HRP)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	HRP
Tested Applications	IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human phospho-GSK-3 Beta (Ser9)
Target	phospho-GSK-3 Beta (Ser9)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Uniprot ID	P49841
Dilution Range	WB=1:500-2000,IHC-P=1:100-500
Expiration Date	12 months from date of receipt.