

Product Datasheet

Tau-4 Antibody (Biotin) (orb1607258)

Description	Tau-4 Antibody (Biotin)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Biotin
Tested Applications	ELISA, IHC-P, WB
Immunogen	Recombinant Human Tau-4 Protein
Target	Tau-4
Preservatives	0.01M PBS(pH7.4)
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Uniprot ID	P10636
Dilution Range	WB=1:500-2000,IHC-P=1:100-500,ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.