

Product Datasheet

human CD86 Antibody (Cy7) (orb1607483)

Description	human CD86 Antibody (Cy7)
Species/Host	Mouse
Reactivity	Human
Conjugation	PE/Cy7
Tested Applications	FC
Target	human CD86
Preservatives	0.01M PBS, 0.5%rAlbumin, 0.03% Proclin300
Form/Appearance	Liquid
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	Mouse IgG1, k
Clonality	Monoclonal
Uniprot ID	035531
Expiration Date	12 months from date of receipt.