

## Product Datasheet

### ZMPSTE24 Immunizing peptide (orb16081)

|                          |                                                                               |
|--------------------------|-------------------------------------------------------------------------------|
| <b>Description</b>       | ZMPSTE24 Immunizing peptide                                                   |
| <b>Preservatives</b>     | 100µg of Lyophilised peptide, formulation<br>De-ionized (DI) water.           |
| <b>Storage</b>           | Shipped at ambient temperature, store at<br>-20°C                             |
| <b>Note</b>              | For research use only                                                         |
| <b>Application notes</b> | This peptide was used as immunogen to<br>produce orb16081 (ZMPSTE24 antibody) |
| <b>Protein Sequence</b>  | ERLQALKTMKQH                                                                  |
| <b>Expiration Date</b>   | 12 months from date of receipt.                                               |