

Product Datasheet

Spectrin alpha chain, non-erythrocytic 1 Antibody (BF647) (orb1610070)

Description	Spectrin alpha chain, non-erythrocytic 1 Antibody (BF647)
Species/Host	Rabbit
Conjugation	BF647
Tested Applications	ICC, IF
Immunogen	KLH conjugated synthetic peptide derived from human Spectrin alpha chain, non-erythrocytic 1
Target	Spectrin alpha chain, non-erythrocytic 1
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q13813
Dilution Range	IF=1:100-500, ICC=1:100-500
Expiration Date	12 months from date of receipt.