

Product Datasheet

Clusterin Antibody (orb1610313)

Description	Clusterin Antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from human CLU
Target	Clusterin
Preservatives	PBS, Glycerol, rAlbumin
Form/Appearance	Lyophilized or Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG/Kappa
Clonality	Recombinant
MW	24/52 kDa
Uniprot ID	P10909
Expiration Date	12 months from date of receipt.