

Product Datasheet

OPA1 Antibody (orb1610321)

Description	OPA1 Antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	IF, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human OPA1
Target	OPA1
Preservatives	1*TBS (pH7.4), 0.05% rAlbumin, 40% Glycerol. Preservative: 0.05% Sodium Azide.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Store at +4°C after thawing. Aliquot store at -20°C. Avoid repeated freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Recombinant
MW	111 kDa
Uniprot ID	O60313
Expiration Date	12 months from date of receipt.