

Product Datasheet

Human Treg Flow Kit (FOXP3 PE/CD4 FITC/CD25 APC) Antibody (orb1610335)

Description	Human Treg Flow Kit (FOXP3 PE/CD4 FITC/CD25 APC) Antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Immunogen	KLH conjugated synthetic peptide derived from human Human Treg Flow Kit (FOXP3 PE/CD4 FITC/CD25 APC)
Target	Human Treg Flow Kit (FOXP3 PE/CD4 FITC/CD25 APC)
Preservatives	0.01M TBS (pH7.4) with 1% rAlbumin, 0.02% Proclin300 and 50% Glycerol.
Storage	Store at 2-8°C.Do not freeze.(12months).
Note	For research use only
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.