

Product Datasheet

LOXL2 Antibody (orb1610336)

Description	LOXL2 Antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human LOXL2
Target	LOXL2
Preservatives	1*TBS (pH7.4), 0.05% rAlbumin, 40% Glycerol. Preservative: 0.05% Sodium Azide.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Recombinant
MW	87 kDa
Uniprot ID	Q9Y4K0
Expiration Date	12 months from date of receipt.