

Product Datasheet

Phospho-PKM2 (Tyr148) Antibody (BF750) (orb1610374)

Description	Phospho-PKM2 (Tyr148) Antibody (BF750)
Species/Host	Rabbit
Reactivity	Human
Conjugation	BF750
Tested Applications	FC, ICC, IF
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human PKM2 around the phosphorylation site of Tyr148
Target	Phospho-PKM2 (Tyr148)
Preservatives	Preservative: 15mM Sodium Azide, Constituents: 1% rAlbumin, 0.01M PBS, pH 7.4
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P14618
Dilution Range	Flow-Cyt=1ug/Test,IF=1:100-500,ICC=1:100
Expiration Date	12 months from date of receipt.