

## Product Datasheet

### Alpha-tubulin Conjugated Antibody (orb1612394)

|                        |                                                               |
|------------------------|---------------------------------------------------------------|
| <b>Description</b>     | Alpha-tubulin Conjugated Antibody                             |
| <b>Species/Host</b>    | Mouse                                                         |
| <b>Conjugation</b>     | Biotin                                                        |
| <b>Storage</b>         | Store at 4°C for 6 months.                                    |
| <b>Note</b>            | For research use only                                         |
| <b>Clonality</b>       | Monoclonal                                                    |
| <b>Source</b>          | Mouse Monoclonal Ab                                           |
| <b>Uniprot ID</b>      | <b>P68366</b>                                                 |
| <b>Dilution Range</b>  | Working with enzyme-conjugated streptavidin, 1: 50 - 1: 1,000 |
| <b>Expiration Date</b> | 12 months from date of receipt.                               |