

Product Datasheet

Beta-tubulin (acetyl-Lys379) Conjugated Antibody (orb1612557)

Description	Beta-tubulin (acetyl-Lys379) Conjugated Antibody
Species/Host	Rabbit
Conjugation	Biotin
Storage	Store at 4°C for 6 months.
Note	For research use only
Clonality	Recombinant
Source	Rabbit Polyclonal Ab
Uniprot ID	P68371, Q13885, Q9BVA1, P04350, P07437, Q13509, Q9BUF5
Dilution Range	Working with enzyme-conjugated streptavidin, 1: 50 - 1: 1,000
Expiration Date	12 months from date of receipt.