

Product Datasheet

MTp53 antibody (FITC) (orb16140)

Description	MTp53 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	FITC
Tested Applications	FC, ICC, IF
Immunogen	KLH conjugated synthetic peptide derived from human p53 around the Mutation site N235K N239Y
Target	Mutant p53 (N235K,N239Y)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	21kDa
Uniprot ID	P04637
Dilution Range	Flow-Cyt=1µg/Test,IF=1:100-500,ICC=1:100
Expiration Date	12 months from date of receipt.