

Product Datasheet

C-type lectin domain family 18 member A Polyclonal Conjugated Antibody (orb1618053)

Description	C-type lectin domain family 18 member A Polyclonal Conjugated Antibody
Species/Host	Rabbit
Conjugation	Biotin
Storage	Store at 4°C for 6 months.
Note	For research use only
Clonality	Recombinant
Source	Rabbit Polyclonal Ab
Uniprot ID	A5D8T8
Dilution Range	Working with enzyme-conjugated streptavidin, 1: 50 - 1: 1,000
Expiration Date	12 months from date of receipt.