

Product Datasheet

MYH-pan Polyclonal Conjugated Antibody (orb1618296)

Description	MYH-pan Polyclonal Conjugated Antibody
Species/Host	Rabbit
Conjugation	Biotin
Storage	Store at 4°C for 6 months.
Note	For research use only
Clonality	Recombinant
Source	Rabbit Polyclonal Ab
Uniprot ID	Q9UKX2 , P13535 , P12883 , Q9Y623 , P13533 , P11055 , P12882
Dilution Range	Working with enzyme-conjugated streptavidin, 1: 50 - 1: 1,000
Expiration Date	12 months from date of receipt.