

## Product Datasheet

# RSK1/2/3 (phospho-Thr573/570/577) Conjugated Antibody (orb1629416)

|                        |                                                               |
|------------------------|---------------------------------------------------------------|
| <b>Description</b>     | RSK1/2/3 (phospho-Thr573/570/577) Conjugated Antibody         |
| <b>Species/Host</b>    | Rabbit                                                        |
| <b>Conjugation</b>     | Biotin                                                        |
| <b>Storage</b>         | Store at 4°C for 6 months.                                    |
| <b>Note</b>            | For research use only                                         |
| <b>Clonality</b>       | Recombinant                                                   |
| <b>Source</b>          | Rabbit Polyclonal Ab                                          |
| <b>Uniprot ID</b>      | <b>P51812, Q15418, Q15349</b>                                 |
| <b>Dilution Range</b>  | Working with enzyme-conjugated streptavidin, 1: 50 - 1: 1,000 |
| <b>Expiration Date</b> | 12 months from date of receipt.                               |