

## Product Datasheet

# RSK1 p90 (Phospho-Ser380) Polyclonal Conjugated Antibody [Out of stock] (orb1629902)

|                        |                                                                               |
|------------------------|-------------------------------------------------------------------------------|
| <b>Description</b>     | RSK1 p90 (Phospho-Ser380)<br>Polyclonal Conjugated Antibody<br>[Out of stock] |
| <b>Species/Host</b>    | Rabbit                                                                        |
| <b>Conjugation</b>     | Biotin                                                                        |
| <b>Storage</b>         | Store at 4°C for 6 months.                                                    |
| <b>Note</b>            | For research use only                                                         |
| <b>Source</b>          | Rabbit Polyclonal Ab                                                          |
| <b>Dilution Range</b>  | Working with enzyme-<br>conjugated streptavidin, 1: 50 -<br>1: 1,000          |
| <b>Expiration Date</b> | 12 months from date of receipt.                                               |