

Product Datasheet

RSK1 p90 (Phospho-Thr359+Ser363) Polyclonal Conjugated Antibody [Out of stock] (orb1629925)

Description	RSK1 p90 (Phospho-Thr359+Ser363) Polyclonal Conjugated Antibody [Out of stock]
Species/Host	Rabbit
Conjugation	Biotin
Storage	Store at 4°C for 6 months.
Note	For research use only
Source	Rabbit Polyclonal Ab
Dilution Range	Working with enzyme-conjugated streptavidin, 1: 50 - 1: 1,000
Expiration Date	12 months from date of receipt.