

Product Datasheet

Trk A (Phospho-Tyr680/681) Polyclonal Conjugated Antibody [Out of stock] (orb1630713)

Description	Trk A (Phospho-Tyr680/681) Polyclonal Conjugated Antibody [Out of stock]
Species/Host	Rabbit
Conjugation	Biotin
Storage	Store at 4°C for 6 months.
Note	For research use only
Source	Rabbit Polyclonal Ab
Dilution Range	Working with enzyme- conjugated streptavidin, 1: 50 - 1: 1,000
Expiration Date	12 months from date of receipt.