

Product Datasheet

Human lambda immunoglobulin light chain, λ -IgLC ELISA Kit (orb1631177)

Description	Human lambda immunoglobulin light chain, λ -IgLC ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.039 μ g/mL-10 μ g/mL. Sensitivity: 0.039 μ g/mL. Assay Time: 1-5 h. Detection wavelength: 450 nm. Sample Volume: 50-100 μ L. Sample Type: serum, plasma.
Expiration Date	Please enquire.